



AGM-4000 Auto A

KEMENKES RI AKL 20101717111



Autocoding



20-600 mg/dL
Rentang pengukuran



0,5 µl/5 detik
Sampel darah sedikit
& hasil cepat

✓ 500 Memory

Strip Auto A

KEMENKES RI AKL 20101718710



Sertifikasi



SPESIFIKASI

Volume sampel	: 0,5 µL
Volume pengukuran	: 5 detik
Kapasitas memori	: 500 memori
Baterai	: 1 baterai Lithium (CR2032), DC 3V
Satuan pengukuran	: mg/dL atau mmol/L
Suhu	: 10 - 40°C, Kelembaban <85%
Ukuran	: 93,5 x 49 x 17,5 (mm)
Layar	: LCD 37 x 42 (mm)
Strip	: Glucodr Auto A Test Strip
Isi	: 25 strip/box dan 50 strip/box
Rentang pengukuran	: 20 - 600 mg/dL (1,1 - 33,3 mmol/L)

Prosedur Tes



Cek Kode Nomor



Masukkan strip
dan sampel darah



Waktu pengukuran
5 detik

Sub-distributor :